


New!

STARfocus

Tubes for Rotor-Gene®



Pack content is sufficient for
1728 reactions



Pack content is sufficient for
3000 reactions



The Rotor-Gene® Style Discs are specially made to be used with the Qiagen Rotor-Gene® Q Real-Time Rotary Analyzer. Two models are available:

- ▶ a 72-well disc with 100 µl tubes (recommended reaction volume 20-25 µl) for use with the original Rotor-Disc® 72 Rotor and Rotor-Disc 72 Locking Ring
- ▶ a 100-well disc with 30 µl tubes (recommended reaction volume 15-20 µl) for use with the original Rotor-Disc® 100 Rotor and Rotor-Disc 100 Locking Ring

The discs are a one-piece "plate" equivalent, having vertically oriented wells compatible with automated reaction setup using a robotic liquid handling system.

Rotor-Gene® Circles (Natural)		
Description	Pack Size	Cat. No.
Disc with 72 x 100 µl Tubes for Rotor-Gene®, recommended reaction volume 20–25 µl	24	I1402-0500
Disc with 100 x 25 µl Tubes for Rotor-Gene®, recommended reaction volume 15–20 µl	30	I1402-0600

STARLAB International GmbH
Neuer Höltigbaum 38 • 22143 Hamburg
Tel.: +49 (0)40 675 99 39 0 • Fax: +49 (0)40 675 99 39 20 • info@starlab.de

Passionate for Science.



www.starlabgroup.com