



UniPix

Innovative laboratory solutions

ATMOS - Pipette Accuracy Tester

Gravimetry

ATMOS

by UniPix



ATMOS is measuring, in dry conditions, the exact displacement volume as well as the leak of the pipette.

The measurement data are stored in the instrument and are shown on the ATMOS, can be printed or can be downloaded on a computer for longitudinal pipette performance monitoring.

ATMOS 1000

ATMOS Multi



CHF 1420, excl VAT

CHF 1790, excl VAT

Distributed in Switzerland by:

Life Systems Design AG

Sagistrasse 8b

8910 Affoltern am Albis

Tel. +41(0)56/6645980

Fax +41(0)56/6645988

mail@lifesystemsdesign.ch

